## **Estimation of Protein by Bradford method**

Dr. Mahesha H B., Yuvaraja's College, Mysore.

Aim: To estimate the protein using Bradford method.

**<u>Principle:</u>** The assay is based on the ability of protein to bind coomassie brilliant blue G250 and form a complex whose *extinction coefficient* is much greater than that of the free dye.

## **Reagents Required:**

**1. Dye Concentrate:** Dissolve 100 mg of coomassie brilliant blue G250 in 50 ml of 95 % ethanol. Add 100 ml of concentrated orthophosphoric acid. Add distilled water to a final volume of 200 ml. store refrigerated in amber bottles; the solution is stable at least 6 months.

- Mix 1 volume of concentrated dye solution with 4 volumes of distilled water for use. Filter with Whatman No. 1 paper.

2. Protein Standard: 100 µg/ml in PBS.

Apparatus and Glass wares required: Test tubes, Pipettes, Colorimeter, etc.,

## **Procedure:**

- 1. Pipette out 0.0, 0.2, 0.4, 0.6, 0.8 and 1 ml of working standard in to the series of labeled test tubes. Also, Pipette out 1 ml of the given sample in another test tube.
- 2. Make up the volume to 1 ml in all the test tubes with PBS. A tube with 1 ml of distilled water serves as the blank.
- 3. Now add 5 ml of diluted dye solution to all the test tubes including the test tubes labeled 'blank' and 'unknown'.
- 4. Mix the contents of the tubes by vortexing / shaking the tubes and allow the colour to develop for at least 5 min but not more than 30 min. The red dye turns blue when it binds protein. Now record the absorbance at 595 nm against blank.
- 5. Then plot the standard curve by taking concentration of protein along X-axis and absorbance at 595 nm along Y-axis.
- 6. Then from this standard curve calculate the concentration of protein in the given sample.

**<u>Result</u>**: The given unknown sample contains ----µg protein/ml.

Volume of	Volume of	Concentration of	Volume of		
standard BSA	distilled water	Protein (µg)	Biuret reagent		
(ml)	(ml)		(ml)	Allow to	A595
0.0	1.0	00	5	develop	0.00
0.2	0.8	1	5	coloue	
0.4	0.6	2	5	from 5	
0.6	0.4	3	5	30 min	
0.8	0.2	4	5		
1.0	0.0	5	5		
1.0 UK	0.0	To be estimated	5		

## **Observations and Calculations**

.