Experiment No 13. Estimation of Succinate Dehydrogenase activity in the eggs/ tissue.

Aim: To estimate succinate dehydrogenase activity in silkworm eggs.

Principle: Succinate dehydrogenase is one of the mitochondrial enzymes, which catalyzes the conversion of succinate to fumerate. In this reaction FAD reduces to FADH₂. In *in vitro*, the lemon yellow colored INT accepts electrons and becomes red coloured farmazan which can be measured at 495 nm.

Reagents Required:

- Sodium phosphate buffer (0.1M, pH 7.4): Mix 16 ml (0.2 M) of monobasic and 84 ml (0.2 M) of dibasic and makeup to 200 ml with distilled water.
- 2. **INT** [2(4-iodophenyl)-3(4-nitrophenyl)-5-Phenyltetrazolium chloride]: 1mg/ml in distilled water.
- 3. Sodium succinate: 15mM
- 4. Glacial acetic acid, Toluene etc.,
- **5. Sample:** Prepare 0.2% egg (5 days old or more) homogenate in cold phosphate buffer/ distilled water using mortar and pestle. Centrifuge the homogenate at 3000 rpm for 10 min, collect the clear supernatant and use as sample.

Procedure for Standard Curve:

- 1. Pipette out 0.0, 0.2, 0.4, 0.6, 0.8 and 1 ml of INT in to the series of labeled test tubes.
- 2. Make up the volume to 1 ml in all the test tubes. A tube with 1 ml of distilled water serves as the blank.
- 3. Now add 1 ml of buffer, 1 ml of sodium succinate and 1 ml of sample (egg homogenate) to all the test tubes including the test tubes labeled 'blank' and 'test'.
- 4. Mix the contents of the tubes by vortexing / shaking the tubes and incubate at 37°C for 24 h or until complete reduction of INT.
- 5. Now add 6 ml of glacial acetic acid to stop the reaction.
- 6. Then add 6 ml of toluene, mix and keep them in a refrigerator for separation of toluene layer.
- 7. Now collect the upper red coloured toluene layer containing farmazan in to a tube, cool to room temperature and record the absorbance at 495 nm against blank.

Reprinted from: A Laboratory Manual on Physiology of Mulberry and Silkworm. Ed. Dr.H.B.Mahesha, Pub. Yuvaraja's College Cooperative Society, University of Mysore, Mysuru, Revised Reprint 2018-19. Then plot the standard curve by taking concentration of farmazan (equivalent to INT) along X-axis and absorbance at 495 nm along Y-axis.

For drawing the standard graph please refer last page figure 1.

Procedure for Test:

- Blank: Take 1 ml of phosphate buffer, 1 ml of sodium succinate, 1 ml of INT and 1 ml of inactivated enzyme sample (distilled water may be used) in a clean dry test tube, and incubate at 37°C for one hour.
- 2. **Test:** Take 1 ml of phosphate buffer, 1 ml of sodium succinate, 1 ml of INT and 1 ml of enzyme sample in a clean dry test tube, and incubate at 37°C for one hour.

After incubation add 6 ml of glacial acetic acid to both blank as well as test to stop the reaction. Then add 6 ml of toluene to each tube, mix and keep them in a refrigerator to separate the red farmazan. After separation, collect the upper red coloured toluene layer containing farmazan in to a cuvette and record the absorbance at 495 nm against blank.

Observations and Calculations:

Standard Curve:

INT (ml)	H ₂ O (ml)	Buffer (ml)	Succinate (ml)	Sample (ml)	Mix, incubate	Acetic acid (ml)	Toluene (ml)	Mix & keep in a	A ₄₉₅
0.0	1.0	1	1	1	at 37°C	6	6	fridge	
0.2	0.8	1	1	1	101 24 11.	6	6	Then	
0.4	0.6	1	1	1		6	6	collect	
0.6	0.4	1	1	1		6	6	upper	
0.8	0.2	1	1	1		6	6	layer	
1.0	0.0	1	1	1		6	6		

Test:

Ex.	INT	Buffer	Succinate	Sample		Acetic	Toluene		A ₄₉₅
	(ml)	(ml)	(ml)	(ml)	Mix &	acid	(ml)	Collect	
					incubate	(ml)		the	
В	1	1	1	1	at 37°C	6	6	farmazan	0.00
Т	1	1	1	1		6	6	as above	

B - Blank

Reprinted from: A Laboratory Manual on Physiology of Mulberry and Silkworm. Ed. Dr.H.B.Mahesha, Pub. Yuvaraja's College Cooperative Society, University of Mysore, Mysuru, Revised Reprint 2018-19. T - Test Optical density of the test: ---SDH activity level = <u>Standard Curve Value X 60</u> Tissue taken (mg) X Incubation time (min)

= ----- μ g of farmazan formed per hour per mg at 37°C.

Report: The SDH activity in the given sample is ----- μg of farmazan formed per hour per mg tissue at 37°C.
